



# **GN7120 Vet**

Akso PCR System

## **GN7120 Vet**

## Akso PCR System





#### Sample In - Result Out

Compact design & All-in-One system Fully automated nucleic acid extraction → PCR amplification → detection → analysis



#### **Quick Test**

Fastest turnaround time of 25 min or less from sample to result



#### **High Throughput**

Dual-module design, 12 test channels Multiple items could be detected simultaneously



#### **Room Temperature Transportation & Storage**

Advanced technology, freeze-dried reagent Convenient and cost-effective transportation



#### **Easy Operation**

Simple workflow Require minimal training and hands-on time



### **Technical Parameters**

Sample Types

Nasopharyngeal swab, conjunctiva swab, anal swab, feces, whole blood, oral and nasal secretion, pleural effusion, ascites

Throughput

12 test channels

Screen

10.4-inch touch screen

**Dimensions** 

420 mm (L)  $\times$  424 mm (W)  $\times$  234 mm (H)

**Storage Capacity** 

10000 data

Working Environment

Temperature: 15~30°C Relative humidity: 10%~70% Air pressure: 70.0~106.0 kPa

Weight

12 kg

#### **Test Items**

∷ 18 Pathogens

Test Items	Pathogen	Specimen	Device
Respiratory Diseases	CDV	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	CIV	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	CPIV	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	CAV- I	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	CAV- II	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	Bb	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	CDV/CPIV/CAV	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	CDV/CPIV/CAV /Bb	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
Digestive Tract	CCV	Feces or Anal Swab	GN7120 Vet
	CPV	Feces or Anal Swab	GN7120 Vet
Diseases	CCV/CPV	Feces or Anal Swab	GN7120 Vet
	Bab	Whole Blood	GN-Pure32/96+GN7120 Vet
Anemia	Bgi	Whole Blood	GN-Pure32/96+GN7120 Vet
	APH	Whole Blood	GN-Pure32/96+GN7120 Vet
	Ehr	Whole Blood	GN-Pure32/96+GN7120 Vet
	Bab/Bgi	Whole Blood	GN-Pure32/96+GN7120 Vet
	APH/Ehr	Whole Blood	GN-Pure32/96+GN7120 Vet
	Bab/Bgi/APH/Ehr	Whole Blood	GN-Pure32/96+GN7120 Vet
Zoonotic Diseases	RV	Oral and Nasal Secretion	GN7120 Vet
	PRV	Oral and Nasal Secretion	GN7120 Vet
	TOX	Feces or Anal Swab	GN7120 Vet
	Lep	Feces or Anal Swab	GN7120 Vet
	Bru	Whole Blood or Serum	GN-Pure32/96+GN7120 Vet
	Нр	Feces or Anal Swab	GN7120 Vet
	TOX/Lep	Feces or Anal Swab	GN7120 Vet
Routine Physical Examination	CDV/CCV/CPV/CAV	Nasopharyngeal and Conjunctiva Swab Feces or Anal Swab	GN7120 Vet

☑: 15 Pathogens

E: Test Items in Development

			. Test items in Developmen
Test Items	Pathogen	Specimen	Device
Respiratory Diseases	FHV	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	FCV	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	MF	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	CF	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	Bb	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	FHV/FCV	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	MF/CF	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	MF/CF/Bb	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	FHV/FCV/MF/CF	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
	FHV/FCV/MF/CF/Bb	Nasopharyngeal and Conjunctiva Swab	GN7120 Vet
Digestive Tract Diseases	FCoV	Feces or Anal Swab	GN7120 Vet
	FPV	Feces or Anal Swab	GN7120 Vet
	TF	Feces or Anal Swab	GN7120 Vet
	GL	Feces or Anal Swab	GN7120 Vet
	GL/TF	Feces or Anal Swab	GN7120 Vet
	FCoV/FPV/TF/GL	Feces or Anal Swab	GN7120 Vet
Feline Infectious	FCoV	Pleural Effusion or Ascites	GN-Pure32/96+GN7120 Vet
	FIPV	Pleural Effusion or Ascites	GN-Pure32/96+GN7120 Vet
Peritonitis	🕵 FCoV/FIPV	Pleural Effusion or Ascites	GN-Pure32/96+GN7120 Vet
Feline Anemia	FeLV	Oral and Nasal Secretion	GN7120 Vet
	FIV	Oral and Nasal Secretion	GN7120 Vet
	FeLV/FIV	Oral and Nasal Secretion	GN7120 Vet
Zoonotic Diseases	TOX	Feces or Anal Swab	GN7120 Vet
	Bh	Feces or Anal Swab	GN7120 Vet
	Di	Feces or Anal Swab	GN7120 Vet
	TOX/Bh	Feces or Anal Swab	GN7120 Vet
Feline Stomatitis	FHV/FCV/FeLV/FIV	Oral and Nasal Secretion	GN7120 Vet



Tel: +86-25-68568508/68568594

Fax: +86-25-68568500

E-mail: sales@getein.com.cn; overseas@getein.com.cn

Web: www.getein.com









